

LEP USB Programming Manual

Purpose:

This document describes the LEP USB interface and how to program it.

Requirements:

- LEP USB Controller
- LEP USB Driver
- USB ready computer running Win98 or Win 2000
- Microsoft Visual C++ version 4.0 or better, or equivalent 32 bit compiler for x86 based system, or Visual basic
- Knowledge of C or Visual basic

Installation:

To install the LEP driver, attach the USB cable to both your computer and the LEP USB controller and power up. The computer should automatically recognize the new hardware. Click the browse button and point to the disk containing the LEP driver software. The driver will be automatically install. You can test the USB interface by running any one of the example programs LEPUSB.exe, Usb_rCfg.exe or Usb_VB.exe. You can find these files on the driver distribution disk or on the LEP web site at www.ludl.com.

Overview:

To communicate with LEP USB Driver one must first enumerate the device. The enumeration of the device returns a device name. This device name is used to open the interface, using CreateFile(). Once you have the handle from CreateFile() you can use ReadFile(), WriteFile(), DeviceIOControl(), and CloseHandle() to communicate to the LEP USB Controller. The hardest part is getting the device name the rest is simply. Example code in console C, Windows Visual Basic, and Windows MFC C++ can be found in appendix B,C, and D respectively.

LEP USB Programming Manual

Device Enumeration:

The device name consists of a number representing a physical port plus the GUID (global unique identifier) appended to it. The GUID for the LEP device is {4d48f140-44e2-11d3-9d64-00e0291dee58}, and a typical device name looks like [\\.\0000000000000002#{4d48f140-44e2-11d3-9d64-0e0291dee58}](#). The device name changes each time you plug in an additional device or plug the device into a different USB port or hub.

There are two ways to get the device name. The easiest method is to read the device name from the registry. When a LEP USB device is plugged in to your computer, the OS detects it and loads the driver. When this driver loads the device name is stored in the registry. The user just has to read the name out of the registry. This method has one disadvantage. It can't be used when more than one USB LEP device is plugged in to your computer, because only the last device name will be recorded in the registry.

To use the registry method simply open the registry key, query the value, and close the registry key. The registry key name is DeviceName and the path is HKEY_LOCAL_MACHINE\System\CurrentControlSet\Services\LepUsb\Parameters\ . Here's an example in Visual Basic.

```
DIM DeviceName as string
DeviceName = GetRegValue(HKEY_LOCAL_MACHINE, _
                        "System\CurrentControlSet\Services\LepUsb\Parameters\", _
                        "DeviceName")

' GetRegValue - Gets the Key value in the registry given a registry key.
Function GetRegValue(hKey As Long, lpszSubKey As String, szKey As String) As String

Dim phkResult As Long, lResult As Long
Dim szBuffer As String, lBuffSize As Long

'Create Buffer
szBuffer = Space(255)          ' Allocate buffer space
lBuffSize = Len(szBuffer)      ' Set the length

RegOpenKeyEx hKey, lpszSubKey, 0, 1, phkResult      'Open the Key, get a handle to it

lResult = RegQueryValueEx(phkResult, szKey, 0, 0, szBuffer, lBuffSize) 'Query the value

RegCloseKey phkResult          'Close the Key

If lResult = ERROR_SUCCESS Then
    GetRegValue = szBuffer      ' return key value
End If
Exit Function
```

The second method to get the device name is to use the device manager. To do this one calls functions in the setupapi DLL. Simply poll the device manager with the LEP GUID for all the devices that match the GUID given. The device manager will return the device names for all the devices currently available. Below is a C example using this method.

Use the DEFINE_GUID macro to build the GUID.

```
DEFINE_GUID(GUID_LEP, 0x4d48f140, 0x11d3, 0x9d, 0x64, 0x00,
            0xe0, 0x29, 0x1d, 0xee, 0x58);
```

LEP USB Programming Manual

This GUID is passed to SetupDiGetClassDevs(), of which returns a handle to the device. The enumeration functions are found in the setupapi library.

```
HDEVINFO hinfo = SetupDiGetClassDevs(&GUID_LEP, NULL, NULL,  
DIGCF_PRESENT | DIGCF_INTERFACEDevice);
```

The first argument identifies the interface you're looking for. The flag bits in the last argument indicate that you are looking for the interfaces exported by the LEP USB device.

Once you have a handle to the device information set, you can perform a enumeration of all the devices that export the particular interface you're interested in. See Microsoft function documentation for more information on setupapi library functions.

Poll the device manager till there are no matching devices left.

```
int i;  
CString Devices[10];          // an array of cstrings  
for (DWORD i=0; ; ++i)  
{  
    SP_INTERFACE_DEVICE_DATA Interface_Info;  
    Interface_Info.cbSize = sizeof(Interface_Info);  
    // Enumerate device  
    if (!SetupDiEnumInterfaceDevice(hInfo, NULL, (LPGUID)  
        &LEP_GUID,i, &Interface_Info))  
    {  
        SetupDiDestroyDeviceInfoList(hInfo);  
        return(i);  
    }  
  
    DWORD needed;              // get the required length  
    SetupDiGetInterfaceDeviceDetail(hInfo, &Interface_Info,  
        NULL, 0, &needed, NULL);  
    PSP_INTERFACE_DEVICE_DETAIL_DATA detail =  
        (PSP_INTERFACE_DEVICE_DETAIL_DATA) malloc(needed);  
    if (!detail)  
    {  
        SetupDiDestroyDeviceInfoList(hInfo);  
        return(i);  
    }  
  
    // fill the device details  
    detail->cbSize = sizeof(SP_INTERFACE_DEVICE_DETAIL_DATA);  
    if (!SetupDiGetInterfaceDeviceDetail(hInfo, &Interface_Info,  
        detail, needed,NULL, NULL))  
    {  
        free((PVOID) detail);  
        SetupDiDestroyDeviceInfoList(hInfo);  
        return(i);  
    }  
  
    char name[MAX_PATH];  
    strncpy(name, detail->DevicePath, sizeof(name));  
    free((PVOID) detail);  
    Devices[i] = name;         // keep a copy of each device name  
} // end of for loop
```

LEP USB Programming Manual

After this code runs you end up with a list of device names, or NULL if no devices could be found ($i = 0$). Each device name will represent one LEP USB controller that is plugged into your computer. If you know that you will only support one LEP USB controller on your system at one time, you can reduce the enumeration code by dropping the for loop and only going through the code once. The device name(s) that are returned from the above code has a port number prefixed to the original GUID. The port number is related to order of the plug and play device on your machine and can not be predetermined. The device name should look like the following.

```
\\.\000000000000000002#{4d48f140-44e2-11d3-9d64-00e0291dee58}
```

This is the complete device name we need to communicate with the LEP USB device.

LEP USB Programming Manual

Device Communications:

Open Device:

To begin communicating with the LEP USB Controller you must first acquire a handle to it. To do this just pass the device name to the CreateFile() function. This is done in the same manner as opening or creating a file. If successful this function will return a handle to the device. If the device is not plugged in, unpowered, or opened by another program this function will fail.

```
HANDLE hUsbDevice = CreateFile(devicename, GENERIC_READ |  
                                GENERIC_WRITE, 0, NULL, OPEN_EXISTING, 0, NULL);
```

Reading and Writing:

Once the handle is acquired you may read and write to the LEP USB device. To read call the ReadFile() function.

```
ReadFile(hUsbDevice, pText, nBytes, &nBytesRead, NULL);
```

The ReadFile() function will try to read the number of byte asked for and return the request in the pointer passed. If the number of bytes requested is larger than what the driver IN buffer has, it will wait till either; the data comes or till it times out. The default timeout is 30 seconds. The driver IN buffer is defaulted to 512 bytes. If this buffer fills up, USB IN packet will be halted till the buffer empties.

To write to the LEP USB device use the WriteFile() function.

```
WriteFile(hUsbDevice, pText, strlen(pText), &nBytesWritten, NULL);
```

The WriteFile() function will write the data passes to the LEP USB device. The function will return when either all the data has been written or an error has occurred. The driver OUT buffer size is created dynamically so any buffer size may be passed.

Close Communication:

When your application has finished using the device, the device should be closed. To do this call CloseHandle() with the device handle. If you do not close the device, you will not be able to access it again without re-setting the device port.

```
CloseHandle( hUsbDevice ) ;
```

LEP USB Programming Manual

Device I/O Controls:

The DeviceIoControl() function allow the application program to read driver status and set driver variable. This function accepts CTL_CODES that are passed to the device driver. Each code represents a function. The following available functions are described below.

Receive Data Length: - IOCTL_LEPUSB_GET_READ_LENGTH

In order to make device reads with out waiting for a time out condition, the application program needs to know how much data exist in the driver's receive buffer. To do this, use the DeviceIoControl function to issue a read IN buffer length command. First construct and fill the ReadLenStruct data structure and call DeviceIoControl() with the data structure as the OUT buffer. When the function returns the data structure is filled with all the IN buffer length data. The data structure returns the following; The length of the currently available data in the buffer. The available free space in the buffer. The length of data up to and including the first line feed (LF) character found. Lastly the length of the data up to and including the first match of the user character passed. If the LF or user character are not found a zero length is returned. To set the user defined character, pass it as the first character in the IN buffer portion of the DeviceIoControl function. See the function definition below.

The DeviceIoControl function sends a control code directly to a specified device driver, causing the corresponding device to perform the specified operation.

```
BOOL DeviceIoControl(
    HANDLE hDevice,           // handle to device of interest
    DWORD dwIoControlCode,    // control code of operation to perform
    LPVOID lpInBuffer,        // pointer to buffer to supply input data
    DWORD nInBufferSize,      // size, in bytes, of input buffer
    LPVOID lpOutBuffer,       // pointer to buffer to receive output data
    DWORD nOutBufferSize,     // size, in bytes, of output buffer
    LPDWORD lpBytesReturned,   // pointer to variable to receive byte count
    LPOVERLAPPED lpOverlapped // pointer to structure for asynchronous operation
);

struct _ReadLenStruct{
    ULONG    DataLength;      // returns length of available data
    ULONG    FreeLength;      // returns the free space in to buffer
    ULONG    LFLength;        // returns the length of data available
                                // till the first Line feed
    ULONG    UserLength;      // returns the length of data available
                                // till the first user defined char.
}    ReadLen;
unsigned long nBytes;
BOOLEAN success;
Unsigned char UserChar = ':'; // set user char to colon ':'

success = DeviceIoControl(hUsbDevice,
    IOCTL_LEPUSB_GET_READ_LENGTH,
    &UserChar, 1, &ReadLen, sizeof(ReadLen),
    &nBytes,
    NULL);
```

LEP USB Programming Manual

Set Timeout: - IOCTL_LEPUSB_TIME_OUT

One can set the driver read time out value by passing an unsigned long value to the time out set function. The units of this value are in milliseconds and the default is 30000 ms.

```
Ulong timeout = 30000; // 30 seconds
success = DeviceIoControl(hUsbDevice,
                          IOCTL_LEPUSB_TIME_OUT,
                          &timeout, 1, &timeout, sizeof(timeout),
                          &nBytes,
                          NULL);
```

Get Driver Version: - IOCTL_LEPUSB_GET_VERSION_INFO

This function will allow an application to read the driver version. This functions returns "LepUsb.sys Version: 1.00.2001.1, Feb 24 2000 16:41:02."

```
void PrintVersion(void )
{
    UCHAR Ver[512];
    ULONG nBytes;
    BOOLEAN Success;

    Success = DeviceIoControl(hUsb, IOCTL_LEPUSB_GET_VERSION_INFO,
                             NULL, 0, Ver, sizeof(Ver), &nBytes, NULL) ;
    if(!Success) printf("ERROR: Read Driver Version.\n");

    Ver[nBytes] = NULL;          // don't forget the NULL terminator
    printf("%s",Ver);           // printf the Version
}
```

Signal Setup: - IOCTL_LEPUSB_SIGNAL_SETUP

This function will allow an application to be signaled when a certain event is satisfied. Events such as; data in receive buffer, line feed or user defined character found, or system not busy can be specified. This function will allow thread program to wait in efficient wait states till the event has occur. Thereby eliminating the constant polling for data. This command is preliminary, contact LEP for more information.

LEP USB Programming Manual

Registry Keys:

The following is a list of registry keys that may be edited. To access the registry run RegEdit.exe. The LEP USB registry keys can be located at:

HKEY_LOCAL_MACHINE\System\CurrentControlSet\Services\LepUsb\Parameters

LEP USB Registry Keys

- DebugLevel Used for debugging, should always be zero.
- InBufferLen The size in bytes of the receive buffer, default 512 bytes.
- TimeOut The time out value in milliseconds, default 30000ms.
- DeviceName This string contains device name of the last LEP USB device loaded. It can be used to open an interface to the LEP USB device.

LEP USB Programming Manual

Appendix A: Related Documents and Web sites:

Universal Serial Bus Specification → www.usb.org

Microsoft Development Network → www.msdn.microsoft.com

Surveying the New Win32 Driver Model for Windows 98 and Windows NT

5.0(Win2000). → www.msdn.microsoft.com/library/periodic/period97/html/wdm1211.htm

Driver Developemnt Kit Website → www.microsoft.com/ddk

LEP USB Programming Manual

Appendix B: Sample code in C:

This example runs in 32-bit console mode. It only scans for one device, then opens it, reads and displays the driver version. Then places the controller in high level mode, writes **ver** and **rconfig** to the LEP controller, wait for the answer, reads the answer, displays the returned string, and closes the handle.

```
// usb_rcfg.cpp : A simple USB DOS console application.
// Author:          Doug Lovett
// Birth:           02/07/2000
// Ludl Electronic Products Copyright 2000

#include "stdafx.h"
#include <objbase.h>
#include <setupapi.h>           // you may have to manually include this
                                library.
#include <initguid.h>
#include <conio.h>

// Defines
#define GUID(LEP_GUID,          // LEP's global, unique identifier
0x4d48f140, 0x44e2, 0x11d3, 0x9d, 0x64, 0x0, 0xe0, 0x29, 0x1d, 0xee, 0x58);

#define CTL_CODE( DeviceType, Function, Method, Access ) ( \
    ((DeviceType) << 16) | ((Access) << 14) | ((Function) << 2) | (Method))

#define METHOD_BUFFERED          0
#define FILE_ANY_ACCESS        0
#define FILE_DEVICE_UNKNOWN    0x00000022

#define LEPUSB_IOCTL_VENDOR_INDEX    0x0800
#define IOCTL_LEPUSB_GET_READ_LENGTH CTL_CODE(FILE_DEVICE_UNKNOWN, \
LEPUSB_IOCTL_VENDOR_INDEX+1,\
                                METHOD_BUFFERED, \
                                FILE_ANY_ACCESS)

#define IOCTL_LEPUSB_GET_VERSION_INFO CTL_CODE(FILE_DEVICE_UNKNOWN, \
LEPUSB_IOCTL_VENDOR_INDEX+2,\
                                METHOD_BUFFERED, \
                                FILE_ANY_ACCESS)

// Function prototypes
int UsbScan(char* );
DWORD GetRxLen(char);
void PrintVersion(void);

// Global Variables
HANDLE hUsb;           // Handle to the USB pipe

void main(int argc, char* argv[])
{
    char    DeviceName[MAX_PATH];
    char    Buffer[1024];
    DWORD   Length, result;
    printf("Usb_rcfg - Version 0.02 %s %s\n", __DATE__, __TIME__);

    // Scan for Usb Device
    if( !UsbScan(DeviceName) ) {
        printf("No devices found!\n");
        return;
    }

    printf("Device Found: %s\n", DeviceName);
}
```

LEP USB Programming Manual

```
hUsb = CreateFile(DeviceName,
                  GENERIC_READ | GENERIC_WRITE,
                  FILE_SHARE_WRITE | FILE_SHARE_READ, NULL,
                  OPEN_EXISTING, 0, NULL);

if( hUsb == INVALID_HANDLE_VALUE ) {
    printf("Error %d: Failed to open USB file handle.\n", GetLastError());
    return;
}
else
    printf("Device Opened: ");

// Print driver version
PrintVersion();

// Put MAC2002 into High Level Mode
sprintf(Buffer, "%c%c", 0xFF, 0x41);
WriteFile(hUsb, Buffer, 2, &Length, NULL);

// Send a 'ver' command
sprintf(Buffer, "ver\n\r");
if(!WriteFile(hUsb, Buffer, strlen(Buffer), &result, NULL))
    printf("ERROR writing to USB.\n");

Length = GetRxLen(':'); // wait for ':'
if(Length) {
    Length = GetRxLen(0); // get total length
    ReadFile(hUsb, Buffer, Length, &result, NULL);
    Buffer[Length] = NULL; // add the NULL terminator
    printf("%s", Buffer);
}

// Send a 'rconfig' command
sprintf(Buffer, "rconfig\n\r");
if(!WriteFile(hUsb, Buffer, strlen(Buffer), &result, NULL))
    printf("ERROR writing to USB.\n");

Length = GetRxLen(':'); // wait for ':'
if(Length) {
    Length = GetRxLen(0); // get total length
    ReadFile(hUsb, Buffer, Length, &result, NULL);
    Buffer[Length] = NULL; // add NULL terminator
    printf("%s", Buffer);
}

CloseHandle(hUsb);
printf("Device Closed.\n");
} // end of main()

// ----- //
// GetRxLen() - If the parameter passed is zero this function simply returns
// the available data in the USB IN buffer. Else this function waits for the
// character given. This work's well for high level commands where a colon(:)
// or linefeed (0x0A) is expected at the end of a incoming transmission.
// For low level mode you should just poll for the length of data that your
// command requires, because in low level command answers are fixed in length.
DWORD GetRxLen(char c)
{
    struct _ReadLenStruct{
        union {
            ULONG    DataLength;
            char    UserChar[sizeof(ULONG)];
        };
        ULONG    FreeLength;
        ULONG    LFLength;
        ULONG    UserLength;
    };
}
```

LEP USB Programming Manual

```
        }
        ReadLen;

ULONG nBytes;
BOOLEAN Success;

while (1) {

    ReadLen.UserChar[0] = c;
    Success = DeviceIoControl(hUsb, IOCTL_LEPUSB_GET_READ_LENGTH,
        &ReadLen, 1, &ReadLen, sizeof(ReadLen),
        &nBytes, NULL) ;
    if(!Success || ( nBytes != sizeof(ReadLen) ) )
    {
        printf("ERROR: Returned from DeviceIoControl.\n");
        CloseHandle(hUsb);
        exit(0);
    }

    if( ReadLen.UserLength || kbhit() ) break;
    if( !c ) break;
}

if( c )        return(ReadLen.UserLength);
else return(ReadLen.DataLength);

}

// -----
// UsbScan() - Scan for first USB device with a GUID that matches the one
//             passed. Return 1 if device found, else 0.
//             This function only detect the first match of the device found.
//             If you have more than one controller attached you will have to
//             scan for the device more than once. See the MFC windows
//             example. You can also get the devicename from the registry.
int UsbScan(char* DeviceName )
{
    // Get handle to the devices
    HDEVINFO hInfo = SetupDiGetClassDevs((LPGUID)&LEP_GUID, NULL, NULL,
        DIGCF_PRESENT | DIGCF_INTERFACEDevice);

    if (hInfo == INVALID_HANDLE_VALUE)
        return(0);

    SP_INTERFACE_DEVICE_DATA Interface_Info;
    Interface_Info.cbSize = sizeof(Interface_Info); // Enumerate device
    if (!SetupDiEnumInterfaceDevice(hInfo, NULL, (LPGUID) &LEP_GUID, 0,
        &Interface_Info))
    {
        SetupDiDestroyDeviceInfoList(hInfo);
        return(0);
    }

    DWORD needed; // get the required length
    SetupDiGetInterfaceDeviceDetail(hInfo, &Interface_Info, NULL, 0,
        &needed, NULL);
    PSP_INTERFACE_DEVICE_DETAIL_DATA detail;
    detail = (PSP_INTERFACE_DEVICE_DETAIL_DATA) malloc(needed);
    if (!detail)
    {
        SetupDiDestroyDeviceInfoList(hInfo);
        return(0);
    }

    // fill the device details
    detail->cbSize = sizeof(SP_INTERFACE_DEVICE_DETAIL_DATA);
    if (!SetupDiGetInterfaceDeviceDetail(hInfo, &Interface_Info, detail,
        needed, NULL, NULL))
    {
        free((PVOID) detail);
    }
}
```

LEP USB Programming Manual

```
SetupDiDestroyDeviceInfoList(hInfo);
return(0);
}

strncpy(DeviceName, detail->DevicePath, MAX_PATH );
free((PVOID) detail);
SetupDiDestroyDeviceInfoList(hInfo);
return(1); // return true, USB device found
}


// ----- //
// PrintVersion() - Gets and printf the driver version information.
void PrintVersion(void )
{
    UCHAR Ver[512];
    ULONG nBytes;
    BOOLEAN Success;

    Success = DeviceIoControl(hUsb, IOCTL_LEPUSB_GET_VERSION_INFO,
        NULL, 0, Ver, sizeof(Ver), &nBytes, NULL) ;
    if(!Success) printf("ERROR: Read Driver Version.\n");

    Ver[nBytes] = NULL; // don't forget the NULL terminator
    printf("%s",Ver); // printf the Version
}

// ----- //
```

Sample Console Output



The screenshot shows an MS-DOS Prompt window with the following text:

```
C:\PROJECTS\dolan\usb_rcfg\Release>usb_rcfg
usb_rcfg - Version 0.02 Feb 22 2000 15:25:34
Device Found: \\.\0000000000000001b#{4d48f140-44e2-11d3-9d64-00e0291dee58}
Device Opened: LepUsb.sys Version: 1.00.2001.1, Feb 16 2000 14:43:41.
Version no. : 4.500 Interface with USB for 2002
:A
Configuration Report

Dev  Address  Label  Id  Description
1      EMOT    X      X  X axis stage  MCMSE
2      EMOT    Y      Y  Y axis stage  MCMSE
3      EMOT    B      B  B aux motor   MCMSE
:A
Device Closed.

C:\PROJECTS\dolan\usb_rcfg\Release>_
```

LEP USB Programming Manual

Appendix C: Sample code in Visual Basic:

This example runs in windows mode. It only scans for one device and opens it. Then the GUI allow the user to enter high level commands. When the program finds an enter key (CR) the program sends the command to the LEP USB device. On a timer event, the length of any available data in the IN buffer is checked. If data is present the data is displayed. This code is very basic and kept to a minimal to illustrate function.

MainForm Code

```
' LEP_USB_VB Example Code
' MainForm 0.001 02/23/2000
' Ludl Electronic Products Copyright 2000 - www.ludl.com
' This code makes function calls to the VB_BAS.bas file.

' Closes the device on exit
Private Sub ExitButton_Click()
Dim Status As Boolean
Status = CloseDevice()      ' close device
End
End Sub

' Initialize the program
Private Sub Form_Load()
Dim Status As Long
Status = OpenDevice()      ' Gets the device name and opens the device
DeviceNameLabel.Caption = VB_USB.lpDeviceName      ' Displays the device Name
Status = WriteDevice((Chr(255) & Chr(65)), 2)      ' Put Controller in High level
mode
If Status = 0 Then          ' Check for errors
    Status = CloseDevice()  ' if errors exit
End
End If
Timer1.Interval = 200      ' Set timer to 200 mSecond interval
Timer1.Enabled = True      ' Enable the timer
End Sub

'Timer1 - Check for data in the IN buffer. If data present, get the data
'and display it in the RxText Box.
Private Sub Timer1_Timer()
On Error GoTo ERROR_HANDLER
Dim Length As Long
Dim ByteRead As Long
Dim RxData As String

Length = GetDeviceLength(0)      ' Get available data in the IN Buffer
If Length <> 0 Then                ' If data exist
    RxData = Space(Length)        ' Allocate the space for the returned string
    BytesRead = ReadDevice(RxData, Length)  ' Get the IN data
    MainForm.RxText.Text = RxData  ' Copy it to the RxText Box
End If
Exit Sub                          ' Exit Sub
ERROR_HANDLER:
MsgBox "RxButton ERROR #" & Str$(Err) & " : " & Error & Chr(13)
End Sub

'KeyPress - Checks to see if the keypressed was an CR. If so, get the command
'and write it to the OUT buffer.
Private Sub TxText_KeyPress(Key As Integer)
On Error GoTo ERROR_HANDLER
Dim TxCommand As String
Dim TxLength As Long

If Key <> 13 Then
    Exit Sub                      ' check for CR, if none exit
End If
```

LEP USB Programming Manual

```
TxCommand = TxText.Text          ' get the command
TxCommand = TxCommand & Chr(13)  ' Add CR for MAC2002 high Level

TxLength = Len(TxCommand)        ' get the length
RxText.Text = TxCommand          ' display the command in the RxText Box

If TxLength <> WriteDevice(TxCommand, TxLength) Then  ' write the text
    MsgBox "TxButton ERROR length not equal"         ' check for error
End If
Exit Sub                                           ' End Sub
ERROR_HANDLER:
    MsgBox "RxButton ERROR #" & Str$(Err) & " : " & Error & Chr(13)
End Sub
```

VB_USB Code

```
' Filename: VB_USB
' Author:   Doug Lovett
' Birth:    02/23/2000
' Ludl Electronic Products Copyright 2000 - www.ludl.com
' Subject:  Functions to access the LEP USB Device.
' Include this file in your VB project and call these functions
' to access the LEP USB device.

' Global data variables
Public hDevice As Long          'Handle to the device
Public lpDeviceName As String 'Copy of the device name

' Registry and File Constants
Const HKEY_LOCAL_MACHINE = &H80000002
Const GENERIC_READ = &H80000000
Const GENERIC_WRITE = &H40000000
Const FILE_SHARE_WRITE = &H2
Const FILE_SHARE_READ = &H1
Const OPEN_EXISTING = &H3
Const CTL_CODE_READ_LEN = &H222004

'Structure Needed for GetDeviceLength
Type GetLenStructure
    DataLength As Long
    FreeLength As Long
    LFLength As Long
    UserLength As Long
End Type

' declare references to external procedures in a dynamic-link library (DLL).
Declare Function RegOpenKeyEx Lib "advapi32" Alias "RegOpenKeyExA" _
    (ByVal hKey As Long, ByVal lpSubKey As String, ByVal ulOptions As Long, _
    ByVal samDesired As Long, phkResult As Long) As Long

Declare Function RegQueryValueEx Lib "advapi32" Alias "RegQueryValueExA" _
    (ByVal hKey As Long, ByVal lpValueName As String, ByVal lpReserved As Long, _
    ByRef lpType As Long, ByVal szData As String, ByRef lpcbData As Long) As Long

Declare Function RegCloseKey Lib "advapi32" (ByVal hKey As Long) As Long

Declare Function CreateFile Lib "kernel32" Alias "CreateFileA" _
    (ByVal lpFileName As String, ByVal dwDesiredAccess As Long, _
    ByVal dwShareMode As Long, ByVal lpSecurityAttributes As Long, _
    ByVal dwCreationDisposition As Long, ByVal dwFlagsAndAttributes As Long, _
    ByVal hTemplateFile As Long) As Long

Declare Function CloseHandle Lib "kernel32" (ByVal hObject As Long) As Boolean

Declare Function DeviceIoControl Lib "kernel32" _
    (ByVal hDevice As Long, ByVal dwIocontrolCode As Long, _
    ByVal lpInBuffer As String, ByVal nInBufferSize As Long, _
    ByRef lpOutBuffer As GetLenStructure, ByVal nOutBufferSize As Long, _
    ByRef lpBytesReturned As Long, ByVal lpOverLapped As Long) As Boolean

Declare Function ReadFile Lib "kernel32" (ByVal hFile As Long, _
    ByVal lpBuffer As String, ByVal nNumberOfBytesToRead As Long, _
    ByRef lpNumberOfBytesRead As Long, ByVal lpOverLapped As Long) As Boolean

Declare Function WriteFile Lib "kernel32" (ByVal hFile As Long, _
    ByVal lpBuffer As String, ByVal nNumberOfBytesToWrite As Long, _
```

LEP USB Programming Manual

```
ByRef lpNumberOfBytesWritten As Long, ByVal lpOverLapped As Long) As Boolean

'OpenDevice - This function reads the device name from the registry and then
'opens the device and stores a handle to the device in hDevice. Returns zero
'on error. This function also stores the full device name in DeviceName.
Function OpenDevice() As Boolean
On Error GoTo ERROR_HANDLER

'Get Device Name from the registry
lpDeviceName = GetRegValue(HKEY_LOCAL_MACHINE, _
                          "System\CurrentControlSet\Services\LepUsb\Parameters\", _
                          "DeviceName", "")
If lpDeviceName = "" Then ' exit on error
    MsgBox "Unable to open device, check connection and power."
    lpDeviceName = "Device Not Found!"
    OpenDevice = 0
    Exit Function
End If
' Try and open the device. This will fail if device not present
hDevice = CreateFile(lpDeviceName, GENERIC_READ Or GENERIC_WRITE, _
                    FILE_SHARE_WRITE Or FILE_SHARE_READ, 0, _
                    OPEN_EXISTING, 0, 0)

If hDevice <= 0 Then ' check for error
    MsgBox "Unable to open device, check connection and power"
    lpDeviceName = "Device Not Found!"
    OpenDevice = 0
End If

Exit Function
ERROR_HANDLER:
    MsgBox "OpenDevice() ERROR #" & Str$(Err) & " : " & Error
End Function

'CloseDevice - Closes the device, always close device after use.
'If you don't close the device after use, you will not be able
'to open it up again without cycle plugging the USB cable.
Function CloseDevice() As Boolean
On Error GoTo ERROR_HANDLER
CloseDevice = CloseHandle(hDevice) ' Close the device
hDevice = 0 ' Null the handle
If CloseDevice = False Then ' Check for errors
    MsgBox "Error closing file" ' Display errors
End If
Exit Function
ERROR_HANDLER:
    MsgBox "CloseDevice() ERROR #" & Str$(Err) & " : " & Error
End Function

' GetDeviceLength - This function returns the data in bytes available
' in the USB IN buffer, if the parameter passed is zero. Else it returns
' the data present up to and including the first match of the parameter
' passed. If there is no match, zero is returned. The later is well suited
' for high level commands were a linefeed of carriage returns is expected
' at the end of an incoming transmission. For low level it would be better
' to poll the data length require, because in low level the data length is
' fixed.
Function GetDeviceLength(c As Byte) As Long
Dim lpDataBytes As Long
Dim nTestByte As Long
Dim TestByte As String
Dim GetLenObject As GetLenStructure

On Error GoTo ERROR_HANDLER
If hDevice <= 0 Then ' check for valid handle
    MsgBox "GetDeviceLength() Handle invalid!"
    Exit Function
End If

If c = 0 Then ' Setup the user search string
    TestByte = ""
    nTestByte = 0
Else
    TestByte = c
    nTestByte = 1
End If

' Call the read length function
If 0 = DeviceIoControl(hDevice, CTL_CODE_READ_LEN, TestByte, nTestByte, _
                      GetLenObject, Len(GetLenObject), lpDataBytes, 0) Then
    MsgBox "GetDeviceLength() Error calling DLL."
```


LEP USB Programming Manual

```
        GetDeviceLength = 0
    Exit Function
End If

If c = 0 Then                                ' Return the Length
    GetDeviceLength = GetLenObject.DataLength
Else
    GetDeviceLength = GetLenObject.UserLength
End If

Exit Function

ERROR_HANDLER:
    MsgBox "GetDeviceLength() ERROR #" & Str$(Err) & " : " & Error
End Function

'ReadDevice - Reads the number of bytes told, from the IN Usb buffer.
'If you ask for more bytes than are available in the IN buffer the
'function will wait till either the data comes or a timeout occurs.
'Returns the amount of data in the returned string.
'Timeout is defaulted to 30sec and may be changed with a DeviceIoControl
'function (see manual) or in the registry at
'(HKEY_LOCAL_MACHINE\System\CurrentControlSet\Services\LepUsb\Parameters\TimeOut).
'The default maximum size for the IN buffer is 512 bytes and can be
'changed with a DeviceIoControl (see manual) or in the registry, at
'(HKEY_LOCAL_MACHINE\System\CurrentControlSet\Services\LepUsb\Parameters\InBufferLen).
Function ReadDevice(lpInBuffer As String, lInLength As Long) As Long
Dim lpBytesRead As Long

On Error GoTo ERROR_HANDLER
If hDevice <= 0 Then                        ' Check for valid handle
    MsgBox "ReadDevice() Handle invalid!"
    Exit Function
End If

If lInLength > Len(lpInBuffer) Then        ' Check string length
    MsgBox "GetDeviceLength() Error in size of string!"
    Exit Function
End If

                                ' Read the Data
If 0 = ReadFile(hDevice, lpInBuffer, lInLength, lpBytesRead, 0) Then
    ReadDevice = 0
Else
    ReadDevice = lpBytesRead                ' Return the length read
End If

Exit Function
ERROR_HANDLER:
    MsgBox "ReadDevice() ERROR #" & Str$(Err) & " : " & Error

End Function

'WriteDevice - Writes the bytes amount told from the string passed to the
'USB OUT Buffer. Returns the amount of bytes written.
Function WriteDevice(lpOutBuffer As String, lOutLength As Long) As Long
Dim lpBytesWritten As Long

On Error GoTo ERROR_HANDLER
If hDevice <= 0 Then                        ' Check for valid handle
    MsgBox "WriteDevice() Handle invalid!"
    Exit Function
End If

                                ' Write the Data to the device
If 0 = WriteFile(hDevice, lpOutBuffer, lOutLength, lpBytesWritten, 0) Then
    WriteDevice = 0
Else
    WriteDevice = lpBytesWritten            ' Return the data amount written
End If

Exit Function
ERROR_HANDLER:
    MsgBox "WriteDevice() ERROR #" & Str$(Err) & " : " & Error
End Function

' GetRegValue - Gets the Key value in the registry given a registry key.
Function GetRegValue(hKey As Long, lpszSubKey As String, szKey As String, _
    szDefault As String) As String
On Error GoTo ERROR_HANDLER
```

LEP USB Programming Manual

```
Dim phkResult As Long, lResult As Long
Dim szBuffer As String, lBuffSize As Long

'Create Buffer
szBuffer = Space(255)          ' Allocate buffer space
lBuffSize = Len(szBuffer)      ' Set the length

'Open the Key
RegOpenKeyEx hKey, lpszSubKey, 0, 1, phkResult

'Query the value
lResult = RegQueryValueEx(phkResult, szKey, 0, 0, szBuffer, lBuffSize)


RegCloseKey phkResult          'Close the Key

'Return obtained value
If lResult = ERROR_SUCCESS Then
    GetRegValue = szBuffer
Else
    GetRegValue = szDefault
End If
Exit Function

ERROR_HANDLER:
MsgBox "GetRegValue() ERROR #" & Str$(Err) & " : " & Error & Chr(13) _
    & "Please exit and try again."

    GetRegValue = szDefault
End Function
```

Sample Visual Basic Output



LEP USB Programming Manual

Appendix D: Sample code in C/C++ MFC Windows:

This example runs in windows mode. It can detect more than one device and mimics a terminal program.

Lep_USB.Cpp - Main window functions

```
// Filename:   UsbLep.Cpp
// Date:      02/23/00
// Author:    Doug Lovett
// Ludl Electronic Products Copyright 2000

// Include files
#include <afxwin.h>
#include <afxcmn.h>
#include <iostream.h>
#include <objbase.h>
#include <initguid.h>
#include "mfc.h"
#include "resource.h"

// The GUID for the LEP USB Device
// {00873FDF-61A8-11d1-AA5E-00C04FB1728B} for LEPUSB.SYS
#define LEP_GUID_STR  "{4D48F140-44E2-11d3-9D64-00E0291DEE58}"

// --- Globals Variables ---
HANDLE hUsb = NULL; // handle to the main usb pipe
CMainEditWin editOb; // Global edit box object
CUsbCode CUsb; // Global Usb function class

// - - - - - //
// Initialize the Application
BOOL CApp::InitInstance()
{
    LPCTSTR cname = AfxRegisterWndClass( 0, NULL, NULL,
                                         LoadIcon(MAKEINTRESOURCE(IDI_ICON1)));

    m_pMainWnd = new CMainWin(cname);
    m_pMainWnd->ShowWindow(m_nCmdShow);
    m_pMainWnd->UpdateWindow();

    // scan for the USD devices
    COpenDeviceDialog diagOb( MAKEINTRESOURCE(IDD_DIALOG1), m_pMainWnd );
    diagOb.DoModal(); // activate modal dialog box
    editOb.SetFocus( );

    return TRUE;
}

CApp App; //Instantiate the application

// - - - - - //
// class CMainWin : public CFrameWnd
// - - - - - //
// The Application's Main Message Map
BEGIN_MESSAGE_MAP( CMainWin, CFrameWnd )
    ON_WM_PAINT()
    ON_WM_SIZE()
    ON_WM_ACTIVATE()
    ON_WM_TIMER()
    ON_COMMAND(IDM_DEVICE_OPEN, OnDeviceOpen)
    ON_COMMAND(IDM_DEVICE_CLOSE, OnDeviceClose)
    ON_COMMAND(IDM_EXIT, OnExit)
    ON_COMMAND(IDM_DEVICE_CLEAR, OnClear)
    ON_COMMAND(IDM_TEST, OnTest)
END_MESSAGE_MAP()

// CMainWin Construrt - Setups up the main window
CMainWin::CMainWin(LPCSTR ClassName)
{
    CString CTitle;
    CTitle = "USB Terminal Program 0.04, ";
    CTitle += (CString)__DATE__ + " " + (CString)__TIME__;
```

LEP USB Programming Manual

```
Create(ClassName, CTitle, WS_OVERLAPPEDWINDOW,
rectDefault, NULL, MAKEINTRESOURCE(IDR_MENU1) );

hUsb = INVALID_HANDLE_VALUE; // Null are USB handle

RECT EditBox; // create an EditBox

structure
    EditBox.left = 0; // set box size
    EditBox.top = 0;
    EditBox.right = GetSystemMetrics(SM_CXSCREEN);
    EditBox.bottom = GetSystemMetrics(SM_CYSCREEN);

// create an edit object
editOb.Create(WS_VSCROLL | WS_HSCROLL | ES_AUTOHSCROLL | ES_AUTOVSCROLL |
ES_LEFT | ES_MULTILINE | ES_NOHIDESEL | ES_WANTRETURN,
    EditBox, this, IDC_EDIT_WIN );
editOb.ShowWindow(SW_SHOW);
editOb.SetFocus( );

CHARFORMAT cf;
editOb.GetDefaultCharFormat( cf );
cf.dwMask = CFM_COLOR;
cf.dwEffects = cf.dwEffects & ~CFE_AUTOCOLOR;
cf.crTextColor = editOb.TextColor = RGB(255,0,0); //RED
editOb.SetWordCharFormat( cf );
editOb.SetDefaultCharFormat(cf);
CTitle += "\n";
editOb.ReplaceSel( CTitle ); // Print the Title

}

// OnTimer Event - On event this function checks for data in the IN USB
// buffer. If data exist, data is read from the device buffer and print
// to the Main Edit Object.
afx_msg void CMainWin::OnTimer( UINT ID )
{
    DWORD length = CUsb.GetLen(); // Get the available data length
    if( length ) // If data? get it
    {
        char *pData = new char[length+1]; // Allocate the buffer
        length = CUsb.Read(pData, length); // Get the data
        pData[length] = NULL; // Don't forget to Null
    }
    terminate it!
    editOb.GreenTextOut(pData); // print it out in green text
    delete(pData); // deallocate the mem
}

// On Main window activation set focus in the edit object
afx_msg void CMainWin::OnActivate(UINT nState, CWnd* pWndOther, BOOL bMinimized)
{
    CWnd::OnActivate( nState, pWndOther, bMinimized );
    if( (nState == WA_ACTIVE) || (nState == WA_CLICKACTIVE) )
    {
        editOb.SetFocus( ); // always place the focus in the edit box
    }
}

// Clear Event - Clears the Edit Objects contents
afx_msg void CMainWin::OnClear()
{
    editOb.SetWindowText("");
}

// OnTest Event - used for debugging to turn off thr timer - Currently unreachable!
afx_msg void CMainWin::OnTest(void)
{
    if( IDYES == MessageBox("Timer?", "Timer", MB_ICONQUESTION | MB_YESNOCANCEL) )
        App.m_pMainWnd->SetTimer(ID_Timer, 30, NULL); //30ms
    else
        App.m_pMainWnd->KillTimer(ID_Timer);
}
```

LEP USB Programming Manual

```
// On Close Event - close device
afx_msg void CMainWin::OnDeviceClose()
{
    CUsb.Close();
}

// Open Device Event - Opens the device and set focus to the main edit object
afx_msg void CMainWin::OnDeviceOpen()
{
    COpenDeviceDialog diagOb( MAKEINTRESOURCE(IDD_DIALOG1), this );
    diagOb.DoModal(); // activate modal dialog box
    editOb.SetFocus( );
}

// Exit event - Send the WM_CLOSE msg that closes the program
afx_msg void CMainWin::OnExit()
{
    int response = IDYES;
    //response = MessageBox("Quit?", "Exit", MB_YESNO );
    if( response == IDYES )
        SendMessage(WM_CLOSE); // terminate program
}

// OnPaint event - Update the screen using contents of virtual window
afx_msg void CMainWin::OnPaint()
{
    CFrameWnd::OnPaint();
}

// On Size event - on size reset the window pos
afx_msg void CMainWin::OnSize( UINT nType, int cx, int cy )
{
    editOb.SetWindowPos( NULL, 0, 0, cx, cy, NULL );
}

// ----- //
// class CMainEditWin : public CEdit
// ----- //
// Main edit box object
BEGIN_MESSAGE_MAP( CMainEditWin, CRichEditCtrl )
    ON_WM_CHAR()
END_MESSAGE_MAP()

CMainEditWin::CMainEditWin()
{
    // make the edit box
    CRichEditCtrl(); // constructor
}

// Onchar event
void CMainEditWin::OnChar(UINT ch, UINT count, UINT flags)
{
    int response = IDYES;
    long nLine;
    char str[1024];

    // Check for CR, if found, send line to USB device
    if(ch == 0x0D)
    {
        nLine = LineFromChar( -1 ); //
        select current line
        response = GetLine( nLine, str, sizeof(str)-2); // get the
        line
        if((str[response-3] == str[response-2]) && (str[response-2] == 0x0D) )
        {
            str[response-2] = 0x0A; // strips off the
            double LF
            str[response-1] = NULL; // Null terminate
        }
        else
            str[response] = NULL; // Null terminate
        CUsb.Write(str); // Send the line to
        the USB device
    }
}
```

LEP USB Programming Manual

```
    }

    editOb.BlueTextOut("");
    CRichEditCtrl::OnChar(ch,count,flags);           // set blue color
                                                    // send char down to the lower class

    // currently edit object is not limited in size, should limit it here
    // if (getnumberlines() > MAXLINES) delete oldest line.

}

// OnCancel event
void CMainEditWin::OnCancel()
{
    KillTimer(ID_Timer);
    CUsb.Close();
    DestroyWindow();
}

// print text to the edit object in red
void CMainEditWin::RedTextOut(char * text)
{
    CHARFORMAT cf;
    GetSelectionCharFormat( cf );
    // check color, change if need be
    TextColor = RGB(255,0,0);           // Red
    if( cf.crTextColor != editOb.TextColor )
    {
        cf.dwMask = CFM_COLOR;
        cf.crTextColor = TextColor;
        SetWordCharFormat( cf );
    }
    editOb.ReplaceSel( text );
}

// print text to the edit object in blue
void CMainEditWin::BlueTextOut(char * text)
{
    CHARFORMAT cf;
    GetSelectionCharFormat( cf );
    // check color, change if need be
    TextColor = RGB(0,0,255);           // Blue
    if( cf.crTextColor != editOb.TextColor )
    {
        cf.dwMask = CFM_COLOR;
        cf.crTextColor = TextColor;
        SetWordCharFormat( cf );
    }
    editOb.ReplaceSel( text );
}

// print text to the edit object in green
void CMainEditWin::GreenTextOut(char * text)
{
    CHARFORMAT cf;
    editOb.GetSelectionCharFormat( cf );
    // check color, change if need be
    editOb.TextColor = RGB(0,100,0);
    if( cf.crTextColor != editOb.TextColor )
    {
        cf.dwMask = CFM_COLOR;
        cf.crTextColor = editOb.TextColor;
        editOb.SetWordCharFormat( cf );
    }
    editOb.ReplaceSel( text );
}

// ----- //
// class COpenDeviceDialog : public CDialog
// ----- //
// Popup window to select USB device

// Open Device Dialog message map
BEGIN_MESSAGE_MAP( COpenDeviceDialog, CDialog )
    ON_COMMAND(ID_OPEN_GUID, OnSelect )
    ON_LBN_DBLCLK( IDC_LIST1, OnSelect )
    ON_COMMAND(IDC_RESCAN, OnReScan )
END_MESSAGE_MAP()
```

LEP USB Programming Manual

```
// Init Dialog box
BOOL COpenDeviceDialog::OnInitDialog()
{
    int i,j=0;
    CString Devices[10];
    CDialog::OnInitDialog();           // call the base class initialize first

    // initialize the edit box
    CEdit* ebpPtr = (CEdit*)GetDlgItem(IDC_EDIT1);           // link the edit box
    Devices[0] = LEP_GUID_STR;                               // get the
device LEP GUID
    ebpPtr->SetWindowText(Devices[0]);           // display the GUID

    i = CUsb.Scan( Devices );           // scan for
the devices
    if(!i)
        // if none found, exit
        {
            MessageBox("Unable to enumerate Usb device","Enumeration Error");
            SendMessage(WM_CLOSE);
            return(0);
        }
    CListBox* lbPtr = (CListBox *) GetDlgItem(IDC_LIST1);           // link the list box
    // init the list box
    while(j<i)
        // fill the list box with
        lbPtr->AddString(Devices[j++]);           // the device
names

    lbPtr->SetCurSel(0);           // default to
the first
    return TRUE;
}

// Open device dialog call backs
// OnSelect button
afx_msg void COpenDeviceDialog::OnSelect()
{
    CString cmsg;
    CListBox* lbPtr = (CListBox*)GetDlgItem(IDC_LIST1);
    int i;
    i = lbPtr->GetCurSel();           // get the index of the selected device name
    if(i == LB_ERR)
        {
            SendMessage(WM_CLOSE); // if error, exit
            return;
        }
    else{
        lbPtr->GetText(i,cmsg);           // get the device name
        CUsb.Open( cmsg );           // open the device
    }
    SendMessage(WM_CLOSE);           // close the dialog box
}

// On Rescan event button - rescan for the device names
afx_msg void COpenDeviceDialog::OnReScan()
{
    int i,j=0;
    CString Devices[10];
    CEdit* ebpPtr = (CEdit*)GetDlgItem(IDC_EDIT1);
    CListBox* lbPtr = (CListBox *) GetDlgItem(IDC_LIST1);

    // Get the GUID
    ebpPtr->GetWindowText(Devices[0]);

    while( lbPtr->DeleteString(0))           // delete all the devices
    i = CUsb.Scan( Devices );
    if(!i)
        {
            MessageBox("Unable to enumerate Usb device","Enumeration Error");
            SendMessage(WM_CLOSE);
        }

    // init the list box
    while(j<i)           // fill the boxes
        lbPtr->AddString(Devices[j++]);
    lbPtr->SetCurSel(0);
}
```

LEP USB Programming Manual

// - - - - - //

UsbCode.Cpp - USB functions

```
//      Filename:      usbcode.cpp
//      Birth:         7/16/99
//      Author:        Doug Lovett
//      Ludl Electronic Products Copyright 2000

// Includes header
#include <afxwin.h>
#include <afxcmn.h>
#include <iostream.h>
#include <afxtempl.h>
#include <objbase.h>
#include <setupapi.h>
#include "devioctl.h"

#include <initguid.h>
#include "resource.h"
#include "mfc.h"

// CTL_CODE defines
#define LEPUSB_IOCTL_INDEX      0x0000
#define LEPUSB_IOCTL_VENDOR_INDEX 0x0800

#define IOCTL_LEPUSB_GET_CONFIG_DESCRIPTOR      CTL_CODE(FILE_DEVICE_UNKNOWN, \
    LEPUSB_IOCTL_INDEX, \
    METHOD_BUFFERED, \
    FILE_ANY_ACCESS)

#define IOCTL_LEPUSB_RESET_DEVICE      CTL_CODE(FILE_DEVICE_UNKNOWN, \
    LEPUSB_IOCTL_INDEX+1, \
    METHOD_BUFFERED, \
    FILE_ANY_ACCESS)

#define IOCTL_LEPUSB_RESET_PIPE      CTL_CODE(FILE_DEVICE_UNKNOWN, \
    LEPUSB_IOCTL_INDEX+2, \
    METHOD_BUFFERED, \
    FILE_ANY_ACCESS)

#define IOCTL_LEPUSB_TIME_OUT      CTL_CODE(FILE_DEVICE_UNKNOWN, \
    LEPUSB_IOCTL_VENDOR_INDEX, \
    METHOD_BUFFERED, \
    FILE_ANY_ACCESS)

#define IOCTL_LEPUSB_GET_READ_LENGTH      CTL_CODE(FILE_DEVICE_UNKNOWN, \
    LEPUSB_IOCTL_VENDOR_INDEX+1, \
    METHOD_BUFFERED, \
    FILE_ANY_ACCESS)

#define IOCTL_LEPUSB_GET_VERSION_INFO      CTL_CODE(FILE_DEVICE_UNKNOWN, \
    LEPUSB_IOCTL_VENDOR_INDEX+2, \
    METHOD_BUFFERED, \
    FILE_ANY_ACCESS)

// LEP USB Device GUID => {4D48F140-44E2-11d3-9D64-00E0291DEE58}
DEFINE_GUID(LEP_GUID,
0x4d48f140, 0x44e2, 0x11d3, 0x9d, 0x64, 0x0, 0xe0, 0x29, 0x1d, 0xee, 0x58);

// external referenced objects
extern CMainEditWin editOb;
extern HANDLE hUsb;
extern CApp App;
```


LEP USB Programming Manual

```
// CUsbCode Object Functions

CUsbCode::CUsbCode()          // Constructor
{
    FirstTimeOpened = 1;    // flag for first time
}

// Scans for all the device that match the GUID passed.
// returns an array of cStrings, containing the device names found,
// and returns the number of device names found. Returns zero if no
// devices found.
int CUsbCode::Scan(CString * Devices)
{
    // Get handle to the devices manager
    HDEVINFO hInfo = SetupDiGetClassDevs((LPGUID)&LEP_GUID, NULL, NULL,
                                         DIGCF_PRESENT | DIGCF_INTERFACEDEVICE);
    if (hInfo == INVALID_HANDLE_VALUE)    // check for error
        return(0);

    for (DWORD i=0; ; ++i)                // loop till all devices are
found
    {
        SP_INTERFACE_DEVICE_DATA Interface_Info;          // declare
        Interface_Info.cbSize = sizeof(Interface_Info);    // Enumerate
        if (!SetupDiEnumInterfaceDevice(hInfo, NULL, (LPGUID) &LEP_GUID, i,
&Interface_Info))
        {
            SetupDiDestroyDeviceInfoList(hInfo);
            return(i);
        }

        DWORD needed;
        // get the required lenght
        SetupDiGetInterfaceDeviceDetail(hInfo, &Interface_Info, NULL, 0, &needed,
NULL);
        PSP_INTERFACE_DEVICE_DETAIL_DATA detail =
(PSP_INTERFACE_DEVICE_DETAIL_DATA) malloc(needed);
        if (!detail)
        {
            SetupDiDestroyDeviceInfoList(hInfo);
            return(i);
        }

        // fill the device details
        detail->cbSize = sizeof(SP_INTERFACE_DEVICE_DETAIL_DATA);
        if (!SetupDiGetInterfaceDeviceDetail(hInfo, &Interface_Info, detail,
needed, NULL, NULL))
        {
            free((PVOID) detail);
            SetupDiDestroyDeviceInfoList(hInfo);
            return(i);
        }

        char name[MAX_PATH];
        strncpy(name, detail->DevicePath, sizeof(name));    // copy the device
name
        free((PVOID) detail);
        // free the mem
        Devices[i] = name;
        // save device name
    }    // end of for loop
}

// Open the device with the given device name
void CUsbCode::Open(CString DeviceName)
{
    // Create the handle to our device
    hUsb = CreateFile(DeviceName,
        GENERIC_READ | GENERIC_WRITE,
        FILE_SHARE_WRITE | FILE_SHARE_READ, NULL,
        OPEN_EXISTING, 0, NULL);

    if( hUsb == INVALID_HANDLE_VALUE )    // process error
    {
        LPVOID lpMsgBuf;
        FormatMessage(
```

LEP USB Programming Manual

```
        FORMAT_MESSAGE_ALLOCATE_BUFFER | FORMAT_MESSAGE_FROM_SYSTEM,
        NULL, GetLastError(),
        MAKELANGID(LANG_NEUTRAL, SUBLANG_DEFAULT), // Default language
        (LPTSTR) &lpMsgBuf, 0, NULL );
    // Display the string.
    MessageBox(NULL,(LPTSTR)lpMsgBuf,"Failed to open Device",MB_ICONHAND);
    LocalFree( lpMsgBuf ); // Free the buffer.
    editOb.SetFocus( );
}
else
{
    if(FirstTimeOpened){                // if this is the first time print
the driver version
        ULONG nBytes;
        if( !DeviceIoControl(hUsb, IOCTL_LEPUSB_GET_VERSION_INFO,
            NULL, 0, msg, sizeof(msg), &nBytes, NULL))
            editOb.RedTextOut( "ERROR: Reading Driver Version.\n");
        msg[nBytes] = NULL;                // don't forget the
NULL terminator
        editOb.RedTextOut( msg);            // print the driver info
        FirstTimeOpened = 0;
    }
    editOb.RedTextOut( "Device opened.\n");
    sprintf(msg,"%c%c",0xFF,0x41);
    Write(msg);                            // switch to
high level (ascii) mode
    CMenu* pCMenu = App.m_pMainWnd->GetMenu();
    pCMenu->EnableMenuItem( IDM_DEVICE_CLOSE,MF_BYCOMMAND | MF_ENABLED);
    pCMenu->EnableMenuItem( IDM_DEVICE_OPEN , MF_BYCOMMAND | MF_GRAYED);
    if(!App.m_pMainWnd->SetTimer(ID_Timer,30,NULL)) //enable the timer for
30ms
        editOb.RedTextOut("Unable to start read length timer!\n");
}

    editOb.SetFocus( );
}

// Close the device,
int CUsbCode::Close()
{
    App.m_pMainWnd->KillTimer(ID_Timer);        // kills the timer
    if( CloseHandle( hUsb) )                    // close the handle to
our device
    {
        editOb.RedTextOut( "Device closed.\n");
        CMenu* pCMenu = App.m_pMainWnd->GetMenu();
        pCMenu->EnableMenuItem( IDM_DEVICE_OPEN, MF_BYCOMMAND | MF_ENABLED);
        pCMenu->EnableMenuItem( IDM_DEVICE_CLOSE, MF_BYCOMMAND | MF_GRAYED);
        App.m_pMainWnd->DrawMenuBar();
        editOb.SetFocus( );
        return(1);
    }
    else
    {
        LPVOID lpMsgBuf;
        FormatMessage(
            FORMAT_MESSAGE_ALLOCATE_BUFFER | FORMAT_MESSAGE_FROM_SYSTEM,
            NULL, GetLastError(),
            MAKELANGID(LANG_NEUTRAL, SUBLANG_DEFAULT), // Default language
            (LPTSTR) &lpMsgBuf, 0, NULL );
        // Display the string.
        MessageBox(NULL,(LPTSTR)lpMsgBuf,"Failed to close Device",MB_ICONHAND);
        LocalFree( lpMsgBuf ); // Free the buffer.
        editOb.SetFocus( );
    }

    return 0;
}

// Write's a string to the usb out port, string must be NULL terminated
// returns number of bytes written
DWORD CUsbCode::Write(char *pText)
{
    DWORD nBytesWritten;

    if(hUsb == INVALID_HANDLE_VALUE)            // check for valid handle
    {
        MessageBox(NULL,"USB OUT Device closed.", "Failed to write to
Device",MB_ICONHAND);
    }
}
```

LEP USB Programming Manual

```
        editOb.SetFocus( );
        return(0);
    }

    if(WriteFile(hUsb, pText, strlen(pText), &nBytesWritten, NULL)) //
write the string
        return nBytesWritten;
    else
    {
        // process errors
        LPVOID lpMsgBuf;
        FormatMessage(
            FORMAT_MESSAGE_ALLOCATE_BUFFER | FORMAT_MESSAGE_FROM_SYSTEM,
            NULL, GetLastError(),
            MAKELANGID(LANG_NEUTRAL, SUBLANG_DEFAULT), // Default language
            (LPTSTR) &lpMsgBuf, 0, NULL );
        // Display the string.
        MessageBox(NULL,(LPTSTR)lpMsgBuf,"Failed to write Device",MB_ICONHAND);
        LocalFree( lpMsgBuf ); // Free the buffer.
        editOb.SetFocus( );
        return(nBytesWritten);
    }
}

// Reads data from the USB IN buffer into the memory pointer passed.
// Memory MUST be allocated for at least the size of the request.
// Returns the amount of data read from the IN buffer.
// If there is not enough data in the IN buffer, the readfile() function
// will wait till either; the data becomes available, or the function times out.
// Default timeout is 30 seconds, and can be adjust with the registry or
// a function call.
DWORD CUsbCode::Read(char * pText, DWORD nBytes)
{
    DWORD nBytesRead;

    if(hUsb == INVALID_HANDLE_VALUE)
    {
        App.m_pMainWnd->KillTimer(ID_Timer);
        MessageBox(NULL,"USB In Device closed.", "Failed to read to
Device",MB_ICONHAND);
        editOb.SetFocus( );
        return(0);
    }

    if(ReadFile(hUsb, pText, nBytes, &nBytesRead, NULL))
        return nBytesRead;
    else
    {
        LPVOID lpMsgBuf;
        FormatMessage(
            FORMAT_MESSAGE_ALLOCATE_BUFFER | FORMAT_MESSAGE_FROM_SYSTEM,
            NULL, GetLastError(),
            MAKELANGID(LANG_NEUTRAL, SUBLANG_DEFAULT), // Default language
            (LPTSTR) &lpMsgBuf, 0, NULL );
        // Display the string.
        MessageBox(NULL,(LPTSTR)lpMsgBuf,"Failed to read Device",MB_ICONHAND);
        LocalFree( lpMsgBuf ); // Free the buffer.
        editOb.SetFocus( );
        return(nBytesRead);
    }

    return 0;
}


// GetLen - This function returns the length of data available in the
// USB IN buffer.
DWORD CUsbCode::GetLen()
{
    struct _ReadLenStruct{
        ULONG    DataLength;
        ULONG    FreeLength;
        ULONG    LFLength;
        ULONG    UserLength;
    }          ReadLen;
    unsigned long nBytes;
    BOOLEAN success;

    success = DeviceIoControl(hUsb, // call the get length function
        IOCTL_LEPUSB_GET_READ_LENGTH,
        NULL, 0, &ReadLen, sizeof(ReadLen),
        &nBytes,
        NULL);
}
```

LEP USB Programming Manual

```
if( nBytes != sizeof(ReadLen) || !success ) // check of errors
{
    App.m_pMainWnd->KillTimer(ID_Timer);
    LPVOID lpMsgBuf;
    FormatMessage(
        FORMAT_MESSAGE_ALLOCATE_BUFFER | FORMAT_MESSAGE_FROM_SYSTEM,
        NULL, GetLastError(),
        MAKELANGID(LANG_NEUTRAL, SUBLANG_DEFAULT), // Default language
        (LPTSTR) &lpMsgBuf, 0, NULL );
    // Display the string.
    MessageBox(NULL,(LPTSTR)lpMsgBuf,"Failed to read device
length!",MB_ICONHAND);
    LocalFree( lpMsgBuf ); // Free the buffer.
    editOb.SetFocus( );
    editOb.RedTextOut("Error reading receive length!\n");
    Close();
    return(0);
}
else
    return(ReadLen.DataLength); // return the available data length
}
// ----- //
```

MFC C/C++ Screen Output example



LEP USB Programming Manual

Appendix E: USB Bandwidth on the MAC2002 USB Controller

The MAC2002 USB controller uses the bulk transfer type USB communication protocol. The communication bandwidth of the MAC2002 USB controller depends on two variables; the first is the data packet size and the second is how much space is available in the USB frame. The latter says that the more devices you have on a USB port, the more the USB frame is used up and the less bandwidth we have. This is obviously an unknown factor, so to simplify things we will assume that the MAC2002 USB Controller is the only device on the USB bus. Just keep in mind that the more USB devices on the bus the less bandwidth we have.

The other factor is the data packet size. The USB protocol sends packet of data out on each USB frame, therefore the larger the data packet the higher the bandwidth. The average data packet size for the MAC2002 USB controller for high level mode is 20 bytes and for low level mode is 6 bytes. So for high level we can expect and bandwidth of 816 Kbytes or better and low level a bandwidth of 352 Kbytes or better. The table below list bandwidth at different data packet sizes.

USB Bulk Transfer Bandwidth

Data Packet Size	Bytes/Second
1	107000
2	200000
4	352000
8	568000
16	816000
32	105000
64	1210000

Note these values are theoretical and only represent the bandwidth of the USB communication protocol in one way. There are other factors that will decrease these values, such as; host computer speed and application load, MAC2002 Controller Interface loading and communication delays, and communication errors (such as in noisy environments).